

TAS2R45-taste receptor, type 2, member 45 Gene

Reference: PTXBC101154

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR59; T2R45; ZG24P; taste receptor type 2 member 45; G-protein coupled receptor 59; taste receptor, type 2, member 45; taste 2 receptor member 45		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 45
GENE ID	259291
ORIGIN SPECIES	Human
ACCESSIONS	BC101154
VECTOR	pDONR223
NCBI SYMBOL	TAS2R45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataactttctgccatcatattttccattctagtggttacatttgtattggaatttgcataatggctcatagcgttggtaattccaccagtggtggaagacaaaagatctcttctgaccaaattgtcactgctctggcggtctccagagttggttctctgggtgtattattaaat
tggtattcaactgtgtgaatccagctttttgtagtgtagaattagaactactgcttataatatctgggcagtaaccggccatttcagcaactggcctgctactagcctcagcatattttattgtcgaagattgccaatcttcccaaccttattttctcgcttaagaggagaggttaagagtgctattct
gggatgctgtgtgggcctttgctatttttggctgtcatctttttgtgtaaacatgaatcagattgtatggacaaaagaatgaaggaacatgacttggaaatcaaaatggcggtgcaatgtacctttcagatacgaactgaaccatgtagcaaaccttagtaccctttactgtaacctgat
cttttctgctgttagtctgttctctgtgtaaacatctcaagaagatgcacctccatggcaaggatcgaagatccagtagccaaggtccacataaaagtgttgcacactgtgatctcttctctgttatgtgcatcttacttgtgtctgtaataatcagtttggagtttgaagaatctggaacaaa
cctgtcttcatgttctgccaagctattggattcagctgttcttcagccaccggttactcctgatttggggaacaagaagctaaagcagacttatcttgcagtttggcgaatgaggtactga