



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TAS2R41-taste receptor, type 2, member 41 Gene

Reference: PTXBC101157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T2R41; T2R59; taste receptor type 2 member 41; taste receptor type 2 member 59; taste receptor, type 2, member 41; taste 2 receptor member 41		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 41
GENE ID	259287
ORIGIN SPECIES	Human
ACCESSIONS	BC101157
VECTOR	pDONR223
NCBI SYMBOL	TAS2R41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]