

C4orf12-chromosome 4 open reading frame 12 Gene

Reference: PTXBC101195

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf12; FBI4; NCRNA00247; WDFY3 antisense RNA 2		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 12
GENE ID	404201
ORIGIN SPECIES	Human
ACCESSIONS	BC101195
VECTOR	pDONR223
NCBI SYMBOL	C4orf12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaaccggtccgcccggccaggcagcggtgctaggcagcaggggcagcgacgcccgcctttccctctctgccccgcctcctcgcgcgcgtcctcgtcgtgggtccccgcgcgcgcgcctcagcctcggcgctcctcgggtttctctccatcaaggccggg
ccgaccccgaggaccatccggaaagtgggggtgtgtgccgtttcccgacgtgtgtggccatctttaatcctcctcctcctcgtgtttccacctcccgctggtgtctgactga