

TMSL1-thymosin-like 1 (pseudogene) Gene

Reference: PTXBC101199

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRWS1; PS1TP5BP1; actin, cytoplasmic 1; PS1TP5-binding protein 1; beta cytoskeletal actin; actin beta		

For research use only.

GENE DESCRIPTION	thymosin-like 1 (pseudogene)
GENE ID	7115
ORIGIN SPECIES	Human
ACCESSIONS	BC101199
VECTOR	pDONR223
NCBI SYMBOL	TMSL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgacaaaccgatatggctgagatggagaaatcgataagtcgaaactgaagaagacagagacgcaagagaaaaatccactgccttccaaagaaacgattgaacaggagaagcaagcaggcgaatcataa