

C15orf37-chromosome 15 open reading frame 37 Gene

Reference: PTXBC101202

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 15 open reading frame 37		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 37
GENE ID	283687
ORIGIN SPECIES	Human
ACCESSIONS	BC101202
VECTOR	pDONR223
NCBI SYMBOL	C15orf37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagaccagcatccgaaagctggtctctatatacaaaagcagaatccctgctacaaagcagtttcccacctccgactctctcccacggcgaacagccctctggtatccaaccaaccaccaacccaatagggccagcaccagagggtccacaatccagagaaggagctacaa
gcccgtgacttgcaagcctcgaaccccaactccaactgacctccctctgcccagagcaagaccacagataccgcccagagcggacctggctggtctacaaaagccaggcgcgccgaccagaggatcctccgggacgtctccgaccggggtgtcccccagccccccc
ggcgcggcgctcaccagctctctccccctctcctag