

AKR1CL1-aldo-keto reductase family 1, member C-like 1 Gene

Reference: PTXBC101204

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	akr; aldo-keto reductase family 1, member C-like 1; aldo-keto reductase family 1 member C1; prostaglandin F synthase		

For research use only.

GENE DESCRIPTION	aldo-keto reductase family 1, member C-like 1
GENE ID	340811
ORIGIN SPECIES	Human
ACCESSIONS	BC101204
VECTOR	pDONR223
NCBI SYMBOL	AKR1CL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgacggatctgaagcaagccattcagtgaggctgaatgatggacccttcacccagtgctgggatttggcacttatgctcctgatcactcccaaaagccaggctgccaggccaccaaaagtggtattgacgtaggcttcgccatattgattcagcactattataccaaaatgaggaggaggttggacaggccatttgggagaagatcgtgatggtaaccgtcaagagagaggaaatattctacaccatcaagctttgggctactttcttcgggcagaattggttcacccggcctagaaggtcactgaagaaactggaccggactatgtagatctcttcattattcatgtaccattgctatgaagggttcttctga