

NANOS3-nanos homolog 3 (Drosophila) Gene

Reference: PTXBC101209

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NANOS1L; NOS3; ZC2HC12C; nanos homolog 3; nanos C2HC-type zinc finger 3		
For research use only.			

GENE DESCRIPTION	nanos homolog 3 (Drosophila)
GENE ID	342977
ORIGIN SPECIES	Human
ACCESSIONS	BC101209
VECTOR	pDONR223
NCBI SYMBOL	NANOS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggacctttgacctgtggacagattacctgggtttggcacacactggttagggctctgagtggaagagggtctgaaaccaggctgagcccccagccagagccagagccaatgctggagccgggtgagccctggagccgatgccagcgcgggagtcggtgccagtccggg
acccaaggatcagaagcgcagcctggagtcctcgcagctccgaacgctgtgctctttctgcaaacacaacggcgagtcctgggcatctaccagtcacgctgctgaaggacgaggctggcagggtgctgtgtccatcctcgggactacgtgtgtccagtcggtccaca
cgtgagcgcgccacacccgacgctctgcccactactggccagggtctacacctcgtctacagccacaccaccgaaactcggcaggcaagaagctgtccggcctgacaaggcgaagacacaggaacacagccaccgccgaggaggaggaggagcaggttcagaggt
gccgggaagtgtgagccttgcctcctgctctccctcatgtccacctag