

ATAD3C-ATPase family, AAA domain containing 3C Gene

Reference: PTXBC101211

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPase family AAA domain-containing protein 3C; ATPase family, AAA domain containing 3A; ATPase family, AAA domain containing 3C		

For research use only.

GENE DESCRIPTION	ATPase family, AAA domain containing 3C
GENE ID	219293
ORIGIN SPECIES	Human
ACCESSIONS	BC101211
VECTOR	pDONR223
NCBI SYMBOL	ATAD3C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggaagatcgacagctgcccgtgtcctggcaggccacggcgatgcctccaaggacggggctctgaccgagggcatgatggacgcctgcgtgcaagactttgtccagcagcaccagcagatgatgcgtggctgaagggggagaggcctggcccgaggacgagcaacc
ctcatcctgagtcctatgggagaccacacccacggagcctggcgcgggacccctcccacccctgcctttgcggcccccgcacatttaggaataactccccgtaataaagtcccacggggccgcaccgctgtgtctattggctgacacggggcggggtttggggcccctaactcccc
ctgggggtcaagggtga