

C17orf54-chromosome 17 open reading frame 54 Gene

Reference: PTXBC101214

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf54; long intergenic non-protein coding RNA 469		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 54
GENE ID	283982
ORIGIN SPECIES	Human
ACCESSIONS	BC101214
VECTOR	pDONR223
NCBI SYMBOL	C17orf54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgtcatcagcaccaccctccagatgtgacaatcagaaatgtgtccagacatccagatgtcccttatggggaggagagaatgccttggtttgcagactccagctggaaccctgctgcaagcagtgaggagcaacagctcggcaatctggaggcaaggtggaggtgga
gaagcacggaacgaagtctctgtctcaaccagcaaaagtgaaccactcccaggtattacctttataatcagacagagcttcaggagccctgtggagtgggcagcatccccagaggccaagggagatattccagacagctctccagcagagaagtftgggaagatgctga