

C12orf36-chromosome 12 open reading frame 36 Gene

Reference: PTXBC101221

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 12 open reading frame 36		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 36
GENE ID	283422
ORIGIN SPECIES	Human
ACCESSIONS	BC101221
VECTOR	pDONR223
NCBI SYMBOL	C12orf36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagattcctaaattcaaaagccagaaggctgggttcctgttccaccctgcctttacctctgtgtgttctctgatgaagacactcatgctccactatttacttacctctgaaacgaagggtgaccagatcagttgttctctgacctgcttgagggaactcagaggctgtggagcaacagatttgggaatgaggaatcgttccctggaagagtcagagcgtctgtagacacattctgctgattgcacgagctcctccctgggaatccactaagactagagaaagaattccttgaggattcagaggctgaagagtgacagacggctttaattgagaaaaagaagcagaaaattgaggagatgtgaggcatcaatgccagcagacgacatgtgacaggtgctag