



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LOC285205-hypothetical LOC285205 Gene

Reference: PTXBC101229

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	hypothetical LOC285205
GENE ID	285205
ORIGIN SPECIES	Human
ACCESSIONS	BC101229
VECTOR	pDONR223
NCBI SYMBOL	LOC285205
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactctgggctctagccaaatgccaaagctatctaaagggaagttgggggttctgtcttggaaaattgcgaatggagccagaaaatcattactgcgcgcggcgggggaagagagagccgccgggagctggtgcgtgacacactgcagctctcgccgggtctccctggactcatcaagcc
cgggtctgcctctggggcgaggggtgagatcataccgcgaacatgcggcccgctctgcgcgcctctccggagagccgcctgtagtactcatccgccattcttaactctgatggaatcttagccctatagagagagggtcaaaaaacaaataaaccaaacacagacaa
acaaaaaataaaccaaatccaattccacccgcacccctccgcacatccgcacgtccacactgcacacttcttaa