

LOC285398-hypothetical locus LOC285398 Gene

Reference: PTXBC101247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical locus LOC285398
GENE ID	285398
ORIGIN SPECIES	Human
ACCESSIONS	BC101247
VECTOR	pDONR223
NCBI SYMBOL	LOC285398
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgtcccatagggcaagtgagcacacacgaagcagtcctcctcccatgtccttctctgtcctggcatccacagaaaggacctcagtcctcatgagctcctctgcagctgctgatcacctcccggaaccctgcacctatcccggtgccacctgacccagact
caagcatgatcctcattctacctgcctcgggaatcacctggattgcaagtcagaacatggctccaggaaccactccaggccttcagggtctgtgtctccaggaggaggctctgagattggctttctggcaagcccacgtggctttgatgccaccaaggttggcatgcccattgtcttct
tcgtga