

C6orf208-chromosome 6 open reading frame 208 Gene

Reference: PTXBC101251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf208; RP1-266L20.3; dJ182D15.1; long intergenic non-protein coding RNA 574		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 208
GENE ID	80069
ORIGIN SPECIES	Human
ACCESSIONS	BC101251
VECTOR	pDONR223
NCBI SYMBOL	C6orf208
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctgtgtcagcctgcacatctccgctcccaccagctgctctcaactgccaggccagactgtggcaaatctcactctcgcgatgtctgggtttccctcgtgtctgggaggtggtgtgcttgggtggcctgagcactgcagtgaatccatgtttccctccagcaccctgtctgtctccaacttggccgacagctctggccaggacgcagcccagctggtgccaccccgacactgttcattcataagacccttgggttctcacttccctcagattttgccaagagaatggtcctgtgtggccagaaggccagcggggtgcagcctgggactgaaagcagagcggctctggtgggaaaaagtag