

RNF212-ring finger protein 212 Gene

Reference: PTXBC101260

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZHP3; ring finger protein 212		
For research use only.			

GENE DESCRIPTION	ring finger protein 212
GENE ID	285498
ORIGIN SPECIES	Human
ACCESSIONS	BC101260
VECTOR	pDONR223
NCBI SYMBOL	RNF212
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaactgggtgttctgtatcgtcttcagccgccccacaggagctgtcttcagctcaccaactcggggcagctgtactgcgacgcctgcctcggcaaggtaaaaagaatgaatgcttgattgttaaagctcctgtctgtacagttttgcttcaagcataccgacgcagatat
ccaggcatttctcatgagcatagacagctgtgtgaagaagtactccagggaacctccagattttagaatttcaagaaaaacacaggaagagattgttagccttctatagagaaagatttctagggttggagaatcccttaggaagtcagtgctgcagatagaacaactacaagtgatgagat
catcacaacaaacagctttcagcacaataaaaagttcagtttcaacaaacacagggatgctgtgtccacctcactcatcagccccgacagactggagtcgatggaagttgatctctccttctccgattagaaaaatctgagatagcagccggcctgcgagaatctccatgattagtc
acctcaagatggacgaatggggccacatctaacagcttcttctgtttcatccatggttgacctgtgttaagccccctgtgccccgagagtggttcatttcagaggttcacatgtttctgcatagacgtctgtcctcactggctgtctccttagcgttcagttctggaaggcaggggaactc
accaactctaaaaacacttccaatatatgctgaggtacaaggccgtctgtttcgttccagcaggctgagggcacctggacacgttcagaacacccgctgtgtctgtgtttctcttccaatttgaagaaaaagtcattttga