

SFXN5-sideroflexin 5 Gene

Reference: PTXBC101312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BBG-TCC; sideroflexin-5; sideroflexin 5		
For research use only.			

GENE DESCRIPTION	sideroflexin 5
GENE ID	94097
ORIGIN SPECIES	Human
ACCESSIONS	BC101312
VECTOR	pDONR223
NCBI SYMBOL	SFXN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgatagacgactacagcatggcgcggtgtagcgctgagcgatgcacgtcttccaactgggcaaacccgcttcagcagacgtctctatggcgcttcaggcactcttgatcatcgaccctcgacactcttgcactgagagacgtctcag
agaggctgtgagctgtggaggactataagcatgggacctgctggcggtgggtcacaatgaacagctctggagtgacagaaatcaagcaggctattctacatccggacacaaatgagaagatcttccatttagaagtcaggttatcttgggacgccaattgtagtcgg
tcttctctggccaaccagacactggcatccactgtctctgagctggctgaaccagagccacaatgcctgtgtaactatgcaaacgcaatgcgacaaagccttcacatccaagttcatccaggatccctggagctgtcatcagcgcgtctccattgctgtggccttaattgcc
tggttcagaaagccaacaagtaccccgagccacccgcttctcatccagaggtttgtgctgttccctgctgtgacccagtgccaatctgcaatgtggtcctgctgacgggtacgggagctggaggaaggagtgatgtctggagcagcgtggcaccctgtggctctcctcaagatcga
gcccagacgacctgctggagacggcgtgacgcgagtggtctgcccagccatcctggtgctacccccgatcgtcatgtccatgctggagaagacggctctcctgagggacgccccggctgctctccctgtgcaaaagcctgtgtgctggcagccttggcctggcctgcc
gctggccatcagcctctcccgcaaatgtcagagattgaacatcccaattagagccggagatagcccagccacgagcagccggacagtggtgtacaacaagggttgtga