

SLC25A35-solute carrier family 25, member 35 Gene

Reference: PTXBC101330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	solute carrier family 25 member 35		
For research use only.			

GENE DESCRIPTION	solute carrier family 25, member 35
GENE ID	399512
ORIGIN SPECIES	Human
ACCESSIONS	BC101330
VECTOR	pDONR223
NCBI SYMBOL	SLC25A35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactcttgatgagtgccctggcagcctgccccctgtgtattaccatccctggaggtgggaagaccagggatgagttgcaaggagaactgcagggccctggcacataccagcggcactaccgaaatgtcttccatgccttcaccatcggaagtgatggccttctg
ccctgcagaaaggcctggccccccctctgtaccagttctgatgaatggcatccgactgggcacatgagctggctgagctggggctacctgcacacagccgaaggcaccacagtcctgccgcagcgcagcagctggggccatggctgggtcatgggagcctacttgg
ggagccccatctacatggtgaagacacacctgcaggcacagcagcctcagaattgctgtagggcaccagtataagcatcagggcatgtttcaggcgcctaaccgagattggccagaacatggtctgttgggttatggcgtgggctctggcggcctgccccgagttatctcgtt
cctccaccagctgtgcaccttctcatccacaggaacctctgagccagtgaggagatcttctccccagagctggaagttggcgtggctgcatgatgagtgccattgcatgttcttggccatggcaccctttgatgtggcctgcacaaggctctacaaccagccacagatgcac
agggcaaggccctcatgtaccgggggatactggacgtctctgtgcagacagctcggaccgagggcattttggcatgtacaaggaaccgagtcaccaagttctcagcaacatcatgcagtgcttcagcgccttactggggaccaagatgttag