

C20orf197-chromosome 20 open reading frame 197 Gene

Reference: PTXBC101349

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf197; chromosome 20 open reading frame 197		
For research use only.			

GENE DESCRIPTION	chromosome 20 open reading frame 197
GENE ID	284756
ORIGIN SPECIES	Human
ACCESSIONS	BC101349
VECTOR	pDONR223
NCBI SYMBOL	C20orf197
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttgccctgtttcagcacagccctattggcaggcagacggctacgggcacagccacaggctgaagtgtcaacacttcaaacacacagacagtacaatgacaaattggagatcagctctaatctcgcccccatttaatgcattgctgaatattcttgaacatagtcacatccacact
gtcccatgacacaagacgctccaaggctgaagatagaggacttctgcagtcaagagagctgggaactcttgacagtcacaatgtcatttgggtattaaaggctctgcaaaagtctgcaccaataaaccccttgattggcttgatccaatgcatgtttccaaaacactacttgcctg
gggacaccttag