



Reference: PTXBC101352

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KV11.1; Kv8.2; RCD3B; potassium voltage-gated channel subfamily V member 2; potassium channel, subfamily V, member 2; potassium channel, voltage gated modifier subfamily V, member 2; voltage-gated potassium channel subunit Kv8.2; potassium voltage-gated c		

GENE DESCRIPTION	potassium channel, subfamily V, member 2
GENE ID	169522
ORIGIN SPECIES	Human
ACCESSIONS	BC101352
VECTOR	pDONR223
NCBI SYMBOL	KCNV2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]