



Reference: PTXBC101363

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALI3; AMIGO-3; amphoterin-induced protein 3; alivin-3; amphoterin-induced gene and ORF 3; adhesion molecule with Ig like domain 3		

For research use only.

GENE DESCRIPTION	adhesion molecule with Ig-like domain 3
GENE ID	386724
ORIGIN SPECIES	Human
ACCESSIONS	BC101363
VECTOR	pDONR223
NCBI SYMBOL	AMIGO3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]