

LOC153684-hypothetical LOC153684 Gene

Reference: PTXBC101385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC153684		

For research use only.

GENE DESCRIPTION	hypothetical LOC153684
GENE ID	153684
ORIGIN SPECIES	Human
ACCESSIONS	BC101385
VECTOR	pDONR223
NCBI SYMBOL	LOC153684
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcacagagctgtgaacaagagaccacggctcgcgaatggcgggtccaggagcctgaaaggggaatgcagcggcgggtgtcaaggacaacattgttttggcgcaaccagcgggtccgtcaccaagaacccgtcgactctgagaaaaagagaagttcggtaccg
agaaactccgtgcagcaagtgtgtgacagcaaaacccggcttcgcgcggcgtggcgaagagccaacggaaacggcgggtgtgtgcgaccgcgaagccggcactagaggcctcgatggagaaagccccgcaccgagacgaggaactgtgcacagcagcactagcagttg
acacaaacagaaaagtgtctcgtctcgtctcgtgggagaaagcgtggatcaaaaccagcaactcagtgagagagaaccccgagccaactga