

H2AFB3-H2A histone family, member B3 Gene

Reference: PTXBC101415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2ABBD; H2AFB; histone H2A-Bbd type 2/3; H2A Barr body-deficient; H2A.Bbd; histone variant H2A, Barr-body deficient; H2A histone family member B3		

For research use only.

GENE DESCRIPTION	H2A histone family, member B3
GENE ID	83740
ORIGIN SPECIES	Human
ACCESSIONS	BC101415
VECTOR	pDONR223
NCBI SYMBOL	H2AFB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaggaggagagacgccgagggtctccggtgctggcggcggggcggaacctgctctgcacccgccgagcggagcttctgtttcagtgagccaggtggagcgcagctacggaggccactacgtcagcgcctgagtcgcacggcgccggtctacctcgctcggt
tattgagtacgtgacggccaaggctcctggagctggcgggcaacgagggccagaacagcggagagcgggaacatcactccctgctgctggacatggtgttcacaacgacaggctactgagcaccctttcaacagcaccacatctcaagtggccctggcgaggactag