

NCRNA00161-non-protein coding RNA 161 Gene

Reference: PTXBC101416

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00161; C21orf100; long intergenic non-protein coding RNA 161		
For research use only.			

GENE DESCRIPTION	non-protein coding RNA 161
GENE ID	118421
ORIGIN SPECIES	Human
ACCESSIONS	BC101416
VECTOR	pDONR223
NCBI SYMBOL	NCRNA00161
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacattattctttgattcttttaacacatgggtgtgggctgaatgggtgtctcccaagaagtcataatgccgggacccctgcccatagtgcctcacacgtgacattgtttggagttagggtgttgaagatgtcattagtcagaatgaggtcacttgagtga