

TMSB4X-thymosin beta 4, X-linked Gene

Reference: PTXBC101790

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTMB4; TB4X; TMSB4; thymosin beta-4; prothymosin beta-4; t beta-4; thymosin, beta 4, X chromosome; thymosin beta 4, X-linked		

For research use only.

GENE DESCRIPTION	thymosin beta 4, X-linked
GENE ID	7114
ORIGIN SPECIES	Human
ACCESSIONS	BC101790
VECTOR	pDONR223
NCBI SYMBOL	TMSB4X
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgacaaaccgatatggctgagatcgagaaatcgataagtcgaaactgaagaagacagagacgcaagagaaaaatccactgcctccaagaaacgattgaacaggagaagcaagcaggcgaatcgtaa