

RNF150-ring finger protein 150 Gene

Reference: PTXBC101992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RING finger protein 150		
For research use only.			

GENE DESCRIPTION	ring finger protein 150
GENE ID	57484
ORIGIN SPECIES	Human
ACCESSIONS	BC101992
VECTOR	pDONR223
NCBI SYMBOL	RNF150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagctcggccacgaccgctggcctgcgacccaacacaaagtcgccccgacccgcgcaagaactgtagccctatccccaggcgaactgcacgtacagggaataagatccggaacgcgttcctgcagaacgcctcagccgtggtcatctcaacgtgggtcc
aacaccaacgagaccatcaccatgccccacgcggtgtagaagacatcgtggccataatgattcctgagccaaaagggaaggagatagtaagcctgctggaagaacatcaccgtgacaatgtacatcaccatcgaacccggaactgacagaatatgtgagccgacttcggtgt
gtttgtcctcatccttcattgtcctgatcatcttccctcgcagtgctcgtcttttattacatccagaggtttcgaatgcaaatgccaggataggaaccagcgcgactgggggatgcagcaaaagaaagccatcagcaaacctcagatcaggaccatcaagaagggtgacaaggaaaca
gagctgtgatttgacaactgtgcagttgtattgaagggtacaagcccaatgacgtgtgctggatcctgcccgcgcgcatctttccacaagtcctgtgttgacccctggctcttagaccatgtacctgtcccatgtgcaagatgaacattcttaagccctagggatccgccaatccgact
gcatggacgacttcccactgacttcgagggctctctggagggtccaccaccaaccagatcacaggtgccagcgacacacagtgaaatgaaagttcagtcactttggacccctgctgctggactgtgggagcctgcaggtgtgccagatagacccatccccaggaggaga
cgctactcttactactaacagtgaattctgccatactaa