

PPP1R11-protein phosphatase 1, regulatory (inhibitor) subunit 11 Gene

Reference: PTXBC102010

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFAP255; HCG-V; HCGV; IPP3; TCTE5; TCTEX5; protein phosphatase 1 regulatory subunit 11; HCG V; hemochromatosis candidate gene V protein; inhibitor-3; protein phosphatase 1, regulatory (inhibitor) subunit 11; protein phosphatase inhibitor 3; t-complex-asso		

For research use only.

GENE DESCRIPTION	protein phosphatase 1, regulatory (inhibitor) subunit 11
GENE ID	6992
ORIGIN SPECIES	Human
ACCESSIONS	BC102010
VECTOR	pDONR223
NCBI SYMBOL	PPP1R11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaggcagggcctgggctgagcgagaccgtcactgagacaacgggtaccgtgacaaccgagcccgagaaacggagccttaccatcaaaactcggaaacggaagccagagaaaaaggtagaatggacaagtgcactgtggacaatgaacacatggggccgctcatcca
aafgctgtgtattatgagaacacctcggcctttggcgagagctccacggaaagtgatgaggaggaagaaggcgtgtgtcatcacactgtgtacgtggccaccgcaaaggacggcgtcgtgcaaccctaggaccgacccccaccacccctccagcctctgacccctccag
ccccctccaggccaatgcagcactaa