

PSMA5-proteasome (prosome, macropain) subunit, alpha type, 5 Gene

Reference: PTXBC102018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PSC5; ZETA; proteasome subunit alpha type-5; epididymis tissue sperm binding protein Li 11n; macropain subunit zeta; macropain zeta chain; multicatalytic endopeptidase complex zeta chain; proteasome (prosome, macropain) subunit, alpha type, 5; proteasome		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, alpha type, 5
GENE ID	5686
ORIGIN SPECIES	Human
ACCESSIONS	BC102018
VECTOR	pDONR223
NCBI SYMBOL	PSMA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttctaccgggtctgagtagcacagggcggtgaatacttttctcccgaagggaagatttcaagtggaaatgccattgaggctatcaagcttggtctacagccattgggatccagacacagagggtgtgtgctagctgtggagaagagaattactccccactgatggagcccgagc
agcattgagaaaaattgtagagattgatgctcacataggttgccatgagtgaggctaatgctgctgaagacctaattgataaagccagagtgaggagacagaaaccactggttcacctacaatgagacaatgacagtgaggagtgtagcccaagctgttccaatctggctttgagttgg
agaagaagatgcagatccagtgccatgtctgtcccttgagtagcattatttggaggagtgatgagaaaggacccagctgtttcatatggacccatctgggaccttggacagtgatgctcagcaattggctctgcttcagaggggtgccagagctcttgcagaagtttaccac
caagctctgactttgaaagaagccatcaagttcactcatctcaacaagtaattggaggagaagctgaatgcaacaacattgagctagccacagtgagcctggccagaattccacatgttcacaaggaagaacttgaagaggttatcaaggacatttaa