

MORN2-MORN repeat containing 2 Gene

Reference: PTXBC102035

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Mopt; MORN repeat-containing protein 2; MORN motif protein in testis; protein containing single MORN motif in testis; MORN repeat containing 2		

For research use only.

GENE DESCRIPTION	MORN repeat containing 2
GENE ID	729967
ORIGIN SPECIES	Human
ACCESSIONS	BC102035
VECTOR	pDONR223
NCBI SYMBOL	MORN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatggttttgaagacttgagcattttcaggagcagtatatgaaggacaatttaaggataatatgtttcatggactggggacttacacattcccaatggggcaagtatactggaaattcaatgaaatagggtggaaggtgaaggggaatatactgatatccaaggactagaatggagt
ggtaactttcattttacagctgctccagacctgaaattaaagcttcacatgtag