

KCNK7-potassium channel, subfamily K, member 7 Gene

Reference: PTXBC103809

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	K2p7.1; TWIK3; potassium channel subfamily K member 7; potassium channel, subfamily K, member 7; potassium channel, two pore domain subfamily K, member 7; two pore domain K ⁺ channel; potassium two pore domain channel subfamily K member 7		

For research use only.

GENE DESCRIPTION	potassium channel, subfamily K, member 7
GENE ID	10089
ORIGIN SPECIES	Human
ACCESSIONS	BC103809
VECTOR	pDONR223
NCBI SYMBOL	KCNK7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggggctaaagccctgggtccgatacgggctcctggtgtggcccaattgctggccctggggcttggggctgtgggtttccaggccctggagggccctcctgcatgcaggctcaggctgagctcagggcagagctggcagccttcaggcagagcataggcctgcctgccac
ccggagctctggaagagctgctgggcactgcccctggccaccaggcccatgggtctccaccctgggcaacagctcagagggcaggacctgggaccttccctcagccctgctcttgcctgccagatcctaccaccacaggttatggccatggcccaatcgcaggcggaaa
ggccttctgcatggtctatgcagccctgggctgcccagcctctagctctgtggccaccctgcgcattgctgctgctgtgctcagccgccacgtgctgggtagcgggtccactggcagctgtcaccggccaggcctgcgctgctgcaggcagttgcactgggactgctggtggc
cagcagcttctgctgctgccagcgtggtgctgtggggccttcagggcgactgcagcctgctggggccgtctacttctgcttcagctcgtcagcaccattggcctggaggacttgcctccggccggccggccagcctgcacccctgattaccactggccagctgcacttctt
ggtggaggacctcactccagggcacggcgtgggagggtag