

COLEC10-collectin sub-family member 10 (C-type lectin) Gene

Reference: PTXBC103816

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CL-34; CLL1; collectin-10; collectin liver protein 1; collectin sub-family member 10 (C-type lectin); collectin-34; collectin subfamily member 10		

For research use only.

GENE DESCRIPTION	collectin sub-family member 10 (C-type lectin)
GENE ID	10584
ORIGIN SPECIES	Human
ACCESSIONS	BC103816
VECTOR	pDONR223
NCBI SYMBOL	COLEC10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatggcttgcaccttgcctcgaagaaccaattatcctcctggtactattctttgcaaatcagagctcgggtcgtgataattgtagccgtcctaccgctgaagtcgtgccacacacacaatttcaccaggacccaaggagatgatggtagaaaggagatccaggagaagaggga
aagcatggcaagtgaggacgcatggggccgaaaggaaattaaaggagaactgggtgatatgggagatcaggcgaatattggcaagactgggccattgggaagaagggtgacaaaggggaaaaagggttgccttggaatactggagaaaaagcgaagcagggtactgtctgtgattgt
ggaagataccggaattgttggaactggaatattgctcggctcaagacatctatgaagttgtcagaatgtgatagcagggaatagggaactgaagagaattctactacatcgtgcaggaagagaactacagggaatccctaaccactgcaggattcgggtggaatg
ctagccatgcccaaggatgaagctgccacacactcgtcgtgactatgttgcgaagtggtcttcttcgggtgttcattggcgtgaatgaccttgaaggaggacagtagatgttcacagacaactccactgcagaactatagcaactggaatgaggggaacccagcgaccct
atggatcagaggactgtgtggagatgctgagctctggcagatggaatgacacagagtgccatctaccatgtactttgtctgtgagttcatcaagaagaaaaagtaa