

MT1P3-metallothionein 1 pseudogene 3 Gene

Reference: PTXBC103840

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	adenylyl cyclase-associated protein 2; 2810452G09Rik; CAP 2; CAP, adenylate cyclase-associated protein, 2 (yeast)		

For research use only.

GENE DESCRIPTION	metallothionein 1 pseudogene 3
GENE ID	140851
ORIGIN SPECIES	Human
ACCESSIONS	BC103840
VECTOR	pDONR223
NCBI SYMBOL	MT1P3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcctgcaccactgggtgctcctgcacgtgcgccgctcctgcaagtgcgaagagtgcaaatgcacctcctgcaagaagagctgctcctgctgcccatgggctgtgcccaagtgtgccagggtgtgtctgcaaggggcgtgcagctgctgtctga