

OSTbeta-organic solute transporter beta Gene

Reference: PTXBC103842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OSTBETA; OSTB; organic solute transporter subunit beta; OST-beta; organic solute transporter beta subunit; solute carrier family 51 subunit beta; solute carrier family 51 beta subunit		

For research use only.

GENE DESCRIPTION	organic solute transporter beta
GENE ID	123264
ORIGIN SPECIES	Human
ACCESSIONS	BC103842
VECTOR	pDONR223
NCBI SYMBOL	OSTbeta
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcacagtgggggctcccggagaccagccgtactgtgtacccagagagctgtggaagagatgcttgggttttctgtgtggaagatgcctccctggaatcattccatccttgcctggcagctgtgtgtcattataagcatgttctcctgggaagaagcatccaggca
agcagaaaagaaaatgcagccaccagaaaagaaactccagaagtctgtcatttggatgaggccaaggatcacaacagcctaacaacctaagagaaacttgcctcagaaaagccaaacttggccaggtggaacttgagttaaagagagagatgtgtcagtttcttcg
gatgtaccagaaactgagagctag