

HMG4L-high-mobility group (nonhistone chromosomal) protein 4-like Gene

Reference: PTXBC103851

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	high-mobility group (nonhistone chromosomal) protein 4-like
GENE ID	128872
ORIGIN SPECIES	Human
ACCESSIONS	BC103851
VECTOR	pDONR223
NCBI SYMBOL	HMG4L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaggcggataagggtgactgtgatagagaatgaaggattatggaccagcctaaggaggcagaacgatcctaagcccccaaaaggccactgtctgattcttctgttctgagaattctgccccagatcaaatccacaaccctggcatctctattggagacgtggcaaa
aaagctgggtgagatgtggaataactaaatgacagtgaaaagcagccttatgtcactaagggtggcaaaagctgaagaagatgagaaggatgtgtgctgactataagtcgaaaggaaagttggacggcacaaaaggtcctgctaaagttgcctgggaaaagatggaagaagaatgaag
aagatggggagggaagagaaggagatgaataa