

## LST1-leukocyte specific transcript 1 Gene

Reference: PTXBC103856

|                |                                                                                                                       |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | B144; D6S49E; LST-1; leukocyte-specific transcript 1 protein; lymphocyte antigen 117; leukocyte specific transcript 1 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | leukocyte specific transcript 1                                                                   |
| <b>GENE ID</b>          | 7940                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC103856                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | LST1                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgttatcgcgaatgatgtaaaggctggagaggagctgggccagggtctcagagcagggaactccatgcatctctgcagaggctgccagtcaggcagtgaggacgtgacctcagggcagagacaagaggcaccgaaggagatccaagagctgactatgcctg  
cattgctgagaacaaaccacctga