



Reference: PTXBC103888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase NHLRC1; EPM2A; EPM2B; MALIN; bA204B7.2; NHL repeat containing 1; NHL repeat-containing protein 1; NHL repeat containing E3 ubiquitin protein ligase 1		

For research use only.

GENE DESCRIPTION	NHL repeat containing 1
GENE ID	378884
ORIGIN SPECIES	Human
ACCESSIONS	BC103888
VECTOR	pDONR223
NCBI SYMBOL	NHLRC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr