



Reference: PTXBC103893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT7.2; HCA2; MAGE-C3; MAGEC4; melanoma-associated antigen C3; MAGE family testis and tumor-specific protein; MAGE-C3 antigen; cancer/testis antigen 7.2; cancer/testis antigen family 7, member 2; hepatocellular carcinoma-associated antigen 2; hepatocellular		

For research use only.

GENE DESCRIPTION	melanoma antigen family C, 3
GENE ID	139081
ORIGIN SPECIES	Human
ACCESSIONS	BC103893
VECTOR	pDONR223
NCBI SYMBOL	MAGEC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]