

IL1F10-interleukin 1 family, member 10 (theta) Gene

Reference: PTXBC103966

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIL1-theta; FKSG75; IL-1HY2; IL-38; IL1-theta; IL1HY2; interleukin-1 family member 10; FIL1 theta; IL-1 theta; IL-1F10 (canonical form IL-1F10a); family of interleukin 1-theta; interleukin-1 HY2; interleukin-1 receptor antagonist FKSG75; interleukin-1 rec		

For research use only.

GENE DESCRIPTION	interleukin 1 family, member 10 (theta)
GENE ID	84639
ORIGIN SPECIES	Human
ACCESSIONS	BC103966
VECTOR	pDONR223
NCBI SYMBOL	IL1F10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgttccctcccatggcaagatactacataaataatgcagaccagaaggctctatcacagaagatggccagctgctggggagatcctgttcagacaaactgctgcagagaagatctgcacacttctaacagagcttgaccgcaccaaggtcccccatttctggggatc
caggagaggagccgctgctggcatgtgtggagacagaagagggccttcctacagctggaggatgtgaacattgaggaaactgtacaagggtggtaagaggccacacgcttcacctcttcagagcagctcaggctccgcttcaggcttgaggctgctgctggcctggctggtt
cctgtgtggccggcagagccccagcagccagctacagctcaccagagagatgagccctcagccgtaccaagttttacttgaacagagctggtag