

IFNA10-interferon, alpha 10 Gene

Reference: PTXBC103972

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFN-alphaC; interferon alpha-10; IFN-alpha-10; interferon alpha-6L; interferon alpha-C; leIF C; interferon alpha 10		

For research use only.

GENE DESCRIPTION	interferon, alpha 10
GENE ID	3446
ORIGIN SPECIES	Human
ACCESSIONS	BC103972
VECTOR	pDONR223
NCBI SYMBOL	IFNA10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtcttttcttactatggcctgtctgtctcagctacaaatccatctgttctctagcgtgtgatctgcctcagaccacagcctgggtaataaggagcccttgatactcctgggacaaatgggaagaatctctcttctcctgcctgaaggacagacatgattccgaatcccca
ggaggagtgtgatggcaaccagttccagaaggctcaagccatctctgtctccatgatgatccagcagaccttcaatctcttcagcagaggactcatctgtctctgggaacagagcctcctagaaaaatttccactgaactttaccagcaactgaatgacctggaagcatgtgtgatac
aggaggttgggtggaagagactccctgatgaatgaggactccatcctgtgtgaggaataacttccaaagaatcactcttataatagaggaaatacagcccttgctcctgggaggtgtcagagcagaatcatgagatccctctcgttttcaacaaacttgcaaaaagattaag
gaggaaggattga