

CMA1-chymase 1, mast cell Gene

Reference: PTXBC103975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYH; MCT1; chymase; alpha-chymase; chymase 1 preproprotein transcript E; chymase 1 preproprotein transcript I; chymase 1, mast cell; chymase, heart; chymase, mast cell; mast cell protease I; chymase 1		

For research use only.

GENE DESCRIPTION	chymase 1, mast cell
GENE ID	1215
ORIGIN SPECIES	Human
ACCESSIONS	BC103975
VECTOR	pDONR223
NCBI SYMBOL	CMA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcttctctctccccctgctgctcttctctgtgctccagagctgaagctgggagatcatcggggcacagaatgcaagccacattccccgccctacatggcctacgtggaattgtaactccaacgggtccctcaaaatttgggtggtttcttataagacggaactttgtgctgac
ggctgctcgttgtgcagggaaggtctataacagtcacccctggagccataacataacagaggaagaagacacatggcagaagcttgaggttataaagcaatccgtcatccaaaataacacttctactctcaccacgatcatgttactaaagtgaaggagaagccagcctgacctg
gctgtggggacactccccctccatcacaattcaactttgcccacctgggagaaatgtgccgggtgctgctggggaagaacaggtgtgtggaagccgggtcagacactctgcaagaggtgaagctgagactcatggatccccagggcctgcagccacttcagagactttgaccacaatc
ttcagctgtgtgtgggcaatcccagggaagacaaaatctgcatttaaggagactctggggccctctctgtgtgctgggggtggccagggcatgctatcctatggacggtcggatgcaaagccccctgctgtcttaccgccgaatctccattaccggccctggatcaaccagatcctgag
gcaaattaa