

C15orf5-chromosome 15 open reading frame 5 Gene

Reference: PTXBC103988

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf5; long intergenic non-protein coding RNA 597		

For research use only.

GENE DESCRIPTION	chromosome 15 open reading frame 5
GENE ID	81698
ORIGIN SPECIES	Human
ACCESSIONS	BC103988
VECTOR	pDONR223
NCBI SYMBOL	C15orf5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaaaaatagggagctacttaaaattttgaacaagcaagtcacatgatgaaatttaattgaggaaacattgtattgtgattgcatatggatggagtgaggacctagaggcaaaagacattagataggagaccaccagaacaaactagacgtgaaagggtgcagctgtggaaatggaaa
gagatgggcaaaacaaaatcaccacagagaaaaaaatcatccagatttcatgatagattggatgtattttgaaacacctctctccatacaaatgaaataa