

C1orf86-chromosome 1 open reading frame 86 Gene

Reference: PTXBC103992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf86; FP7162; Fanconi anemia core complex-associated protein 20; FANCA-associated protein of 20 kDa; Fanconi anemia-associated protein of 20 kDa; Fanconi anemia core complex associated protein 20		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 86
GENE ID	199990
ORIGIN SPECIES	Human
ACCESSIONS	BC103992
VECTOR	pDONR223
NCBI SYMBOL	C1orf86
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcggcgccggagccgctgggtgagccgagggccgccccggcgggcgggccttctggcgccgccccgtgttctcctgggggtgatgagcgggagcggctctgggcccagctactgcgcacgggtgagccccgagctgatcctggatcacgaggtgcc
ttcactgcccgccttcccaggacaggagccccaggtgcccggccggagcccaactgaagtcttcaactgtcggaccacaagaccttttctggacacctttcccgccgacctgtggggcccgggccgttctaccggctgcttcaaggggcaggaggccacctggaatccccgccaggtcc
ctgccccagcggcccgacctgatccctgcagggccccaggggtggagcagcagccgtctgtggagggtgccgcccctgcgcagctgccccatgtgccagaaggagttcgccccaggtgaccagctggatgttgacagccacctggccagtgcttggccgaaagcacag
aagacgtgacgtgtcttcatggatgtgcaagacttgaatgacctggcagtcagggtga