

YIPF7-Yip1 domain family, member 7 Gene

Reference: PTXBC103996

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein YIPF7; FinGER9; YIP1 family member 7; five-pass transmembrane protein localizing in the Golgi apparatus and the endoplasmic reticulum 9; Yip1 domain family member 7		

For research use only.

GENE DESCRIPTION	Yip1 domain family, member 7
GENE ID	285525
ORIGIN SPECIES	Human
ACCESSIONS	BC103996
VECTOR	pDONR223
NCBI SYMBOL	YIPF7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatttgcctaaattctcacacaaagctacatttactggagatctttctataaaaaataacagaggatgtcaaaactggcacaatttgactctgattttaccaatctaattttactattgataaccaggagcagagtggtaatgactctaagcctatggaatctttatggatctagaaaacaa
caagctggtgagcagcctcagcctgcctctttgttccatcagagatgctcatgtcatcggttacgcaggacaatttttcagccagcatccaactcagattattattcacaatctcttaccattgacagttttgatgaagagcctcctttgctagaagataagaactggaatccattttgatcacata
tggcaaaaaacttga