

PRSS2-protease, serine, 2 (trypsin 2) Gene

Reference: PTXBC103997

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRY2; TRY8; TRYP2; trypsin-2; anionic trypsinogen; protease serine 2 preproprotein; protease, serine, 2 (trypsin 2); trypsin II; trypsinogen 2; protease, serine 2		

For research use only.

GENE DESCRIPTION	protease, serine, 2 (trypsin 2)
GENE ID	5645
ORIGIN SPECIES	Human
ACCESSIONS	BC103997
VECTOR	pDONR223
NCBI SYMBOL	PRSS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccactcctgatccttctgttgcagctgctgttgcctcccttggatgatgacaagatcgttggggctacatctgtgaggagaattctgtccctaccaggtgtccttgaattctggctaccacttctgcgggtgctccctcagcgaacagtggtggtgtcagcaggtca
ctgctacaagtcctccatccaggtgagactgggagagcacaacatcgaagtcctggagggaatgaacagttcatcaatgcagccaagatcatccgccaccccaatacaacagccggactctggacaatgacatcctgctgataagctctcctcactgccgtcatcaattcccgcgtg
tccgccatctctgcccactgcccctccagctgctggcaccgagtcctcctcctcctggctggggcaacactctgagttctgtgcccactaccagacgagctgcagtgccctggatgctcctgctgagccaggtgagtggaagcctcctaccctggaagattaccaacaacatgttc
tgtgtgggtctcctcagggaggcaaggattcctgccaggggtgattctgtgtgcccctgtggtctccaatggagagctccaaggaaattgtctcctggggctatggctgtgccagaagaacagcgctggagtctaccaaggtctacaactatgtggactggattaaggacacatagctgc
caacagctaa