

C6orf195-chromosome 6 open reading frame 195 Gene

Reference: PTXBC104008

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf195; bA145H9.2; chromosome 6 open reading frame 195		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 195
GENE ID	154386
ORIGIN SPECIES	Human
ACCESSIONS	BC104008
VECTOR	pDONR223
NCBI SYMBOL	C6orf195
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttctatccacttgaccttcaggacattccctggaagcaggggaagtgttgccttcctgtctccgagggagaaaggcatttgatggcatggagccactctgccagcaggagccagggccgcttaagggggctacagcatgcctcagattgctacaggctcctctcccctccc
gggtccggccttgctggcaccacccctcagttcctgccccatctcccactgcgggtgtctgacagcctggagcatttctccttcactttcacgggtccacacccctcaagatcaacagtgaccaggccactccctgccatctgagaagcccagaaaccatctccgtcagaggagatt
agcaatgaagcttctga