

C1orf157-chromosome 1 open reading frame 157 Gene

Reference: PTXBC104011

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf157; NCRNA00303; long intergenic non-protein coding RNA 303		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 157
GENE ID	284573
ORIGIN SPECIES	Human
ACCESSIONS	BC104011
VECTOR	pDONR223
NCBI SYMBOL	C1orf157
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtataagtcttggacttgggagacccaaaagagaagggaaggtgttggcaatctctgaagagactaaagaactcagaacaactgaacaagccagcagagccatcaactcacccctacagtcctacatacacagactccatgaaagcattggccatcagctcactggtttctgcc
ccaaattcacacactccctctatcacaaagccaccgcccacgctggcctctcgaatgctgtggaggagagagggcagacaacgttgcctccctgactcccatgtgcgggttcaaacactcagtgga