



Reference: PTXBC104037

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSNB1F; FIGLER4; leucine-rich repeat, immunoglobulin-like domain and transmembrane domain-containing protein 3; fibronectin type III, immunoglobulin and leucine rich repeat domains 4; leucine-rich repeat, immunoglobulin-like and transmembrane domains 3; 1		

For research use only.

GENE DESCRIPTION	leucine-rich repeat, immunoglobulin-like and transmembrane domains 3
GENE ID	345193
ORIGIN SPECIES	Human
ACCESSIONS	BC104037
VECTOR	pDONR223
NCBI SYMBOL	LRIT3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]