

BRP44-brain protein 44 Gene

Reference: PTXBC104157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRP44; mitochondrial pyruvate carrier 2; brain protein 44		
For research use only.			

GENE DESCRIPTION	brain protein 44
GENE ID	25874
ORIGIN SPECIES	Human
ACCESSIONS	BC104157
VECTOR	pDONR223
NCBI SYMBOL	BRP44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgccgcgggtgccgaggcctgcgggccacctaccaccggctctcgataaagtggagctgatgtgcccagaaattgagggcgttgtacaaccatccagcagggtccagaacagttttctctgggtccaattatgaatgggggttggtgtgtgctggattggctgatatggc
cagacctgcagaaaaacttagcacagctcaatctgctgtttgatggctacagggtttatttggtcaagatactcactgttaattattccaaaaaattggagctgtttgctgttaattcttttgggggcagcaggagcctctcagcttttctgtattggagataaaccaagaactaaaagctaaag
cacacaataa