



Antibodies and Proteins

# ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

## C3orf36-chromosome 3 open reading frame 36 Gene

Reference: PTXBC104162

|                               |                                                                     |                     |                |
|-------------------------------|---------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | uncharacterized protein C3orf36; chromosome 3 open reading frame 36 |                     |                |
| <b>For research use only.</b> |                                                                     |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 3 open reading frame 36                                                               |
| <b>GENE ID</b>          | 80111                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC104162                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | C3orf36                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]