

LOC441493-hypothetical gene supported by BC067869 Gene

Reference: PTXBC104165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical gene supported by BC067869
GENE ID	441493
ORIGIN SPECIES	Human
ACCESSIONS	BC104165
VECTOR	pDONR223
NCBI SYMBOL	LOC441493
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagttacggaagtgggtcttgggacgttcgtaccctcggcggcgtgtcaggtgggaacatggccgcgccacgtcgcgtctcccatgggatacaccatgttgaggagtgccatccgtattgaggatgaggacacctgaactcacagagcaggctgtgagttta
gagctgtctctaaactcaggtggagctgtccctgtgtggcagggagcaggggtgggacgtttctctggcggcgttaactggaggaggggtctttggacctgtagtctatggatctgtgtgataccgctccgagaaactgaaggccctcttggaaagaatcagggcgccgcg
gctatctcctggcagataccacaggatgagaagcggggagatgtcaggctgtttacaggttctctagggacataccacaagggaggaggaggaggaatcatgtgtatgcaagccctcctaggatagccatgaggtgggaaaaaggagaactcagagtgctgggtccctc
ctgtgtgctgtggggcttggatgcgaggtttgaatgtcagcaatttaccttctccactccatctacccagctaaaaagggactacgtagtaccataactaataatttttaatgtgtccttagtccaatgcaatgtcaataaaaaatcttaggaactttatttttttctcatgaagta
a