



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MC2R-melanocortin 2 receptor (adrenocorticotrophic hormone) Gene

Reference: PTXBC104169

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACTHR; adrenocorticotrophic hormone receptor; ACTH receptor; MC2 receptor; adrenocorticotropin receptor; corticotropin receptor; melanocortin 2 receptor (adrenocorticotrophic hormone); melanocortin 2 receptor		

For research use only.

GENE DESCRIPTION	melanocortin 2 receptor (adrenocorticotrophic hormone)
GENE ID	4158
ORIGIN SPECIES	Human
ACCESSIONS	BC104169
VECTOR	pDONR223
NCBI SYMBOL	MC2R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacacatattcaactctgtatgaaacatcaacaacacagcaagaataattccgactgtctctgtgtgttttccggagagagatattttcaacaattccattgttggaatttggagaactgatctctctgtgctgtgttcaagaataaagaattccaggcaccatgtacttttcatctgtacttggccatctatgatctgtggcgacctatataaagacttggaaaataatctgatcatattgagaacaatggcgctatctcaaggccactgtggcagtttgaaacacacagccgatgacatcatgacatccctgtttgtctctccctctgttggctccactctcagcctgtctgtgattgtctggga
ccgtacatcacatcattctccagcactctcggtacacacagcatctgacatctgcgcacatctgtgtgtgtgtcttgcactgtcaltctggacgtttctgcagcgggactggcatcacatgtgtgatcttcccatcatgtggccacagtgatcaccttcaactgtcgttctccgtgatctgtgttctcat
cctgtgctctatgtgcacatgtctctgtggcgtgacatccacacagggaagatctcaccctccccagagccaacatgaaggggccacacatgacacatctcaggtgtgttcttctctgtctgtgcttcatctgtgatcatattctcccaagtaaccttactgc
gcctctacatctctctctccaggtgacccatctgtatcatctgtcaatcgctcattgaccctctatatactctccggaagccagagctcaggcagcattcaaaaagatgaltcttcagcaggtactctgtatg

More info: Please contact us for any question: contact@proteogenix.fr