

LOC255025-hypothetical LOC255025 Gene

Reference: PTXBC104171

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC255025
GENE ID	255025
ORIGIN SPECIES	Human
ACCESSIONS	BC104171
VECTOR	pDONR223
NCBI SYMBOL	LOC255025
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtactggaaaacttctaataatcctgtcatltaagaagaagtagtattatcagaagaatcgtagaagtagaactgtgtgtgtcagtttaacaaaagaagaacacattgctttccaactggcccccggccattttccagttgaacatcaaatattttgtgtcttcaactgccaataaa
gggcaacagaaatggtagaggtccctctgcagtggaattacccttatgtctttctgtctatgctgatccacctgacttgaactgggtctattccatggggaattggaacatggggcaattaacactttctccactttatcccatcaacaggaaatgaaatataaccaactaactcagagagac
atcatggaagatgcacacagtgaatga